


ABD-1970

## Chemical Properties

CAS No. : 2010154-82-0  
 Formula: C<sub>21</sub>H<sub>24</sub>ClF<sub>6</sub>N<sub>3</sub>O<sub>3</sub>  
 Molecular Weight: 515.88  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



## Biological Description

Description	ABD-1970 is a highly potent and selective MGLL inhibitor.
Targets(IC50)	Others, MAGL

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9384 mL	9.6922 mL	19.3844 mL
5 mM	0.3877 mL	1.9384 mL	3.8769 mL
10 mM	0.1938 mL	0.9692 mL	1.9384 mL
50 mM	0.0388 mL	0.1938 mL	0.3877 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481