


DL-Arginine

Chemical Properties

CAS No. :	7200-25-1
Formula:	C ₆ H ₁₄ N ₄ O ₂
Molecular Weight:	174.2
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	DL-Arginine is deployed in the study of amino acid complexation dynamics and the formation of crystal structures.
Targets(IC50)	Others,Endogenous Metabolite

Solubility Information

Solubility	H ₂ O: 10 mM,Sonication is recommended. DMSO: Insoluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.7405 mL	28.7026 mL	57.4053 mL
5 mM	1.1481 mL	5.7405 mL	11.4811 mL
10 mM	0.5741 mL	2.8703 mL	5.7405 mL
50 mM	0.1148 mL	0.5741 mL	1.1481 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481