


Brexpiprazole

Chemical Properties

CAS No. :	913611-97-9
Formula:	C ₂₅ H ₂₇ N ₃ O ₂ S
Molecular Weight:	433.57
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Brexpiprazole (OPC-34712) is a partial agonist of human 5-hydroxytryptamine (5-HT) 5-HT _{1A} and dopamine D ₂ receptors.
Targets(IC ₅₀)	5-HT Receptor, Adrenergic Receptor, Dopamine Receptor
In vitro	Brexpiprazole is a potent partial agonist at human 5-hydroxytryptamine (5-HT) 5-HT _{1A} (K _i =0.12 nM) and dopamine D _{2L} (K _i =0.3 nM) receptors, and an antagonist at 5-HT _{2A} receptors (K _i =0.47 nM). It also shows potent antagonist activity at human noradrenergic α _{1B} (K _i =0.17 nM) and α _{2C} receptors (K _i =0.59 nM). Furthermore, this drug displays moderate affinity for human D ₃ , 5-HT _{2B} and 5-HT ₇ receptors, as well as α _{1A} , and α _{1D} adrenergic receptors. Brexpiprazole potentiated NGF-induced neurite outgrowth in PC12 cells. It could significantly potentiate the effects of fluoxetine (or paroxetine) on neurite outgrowth[1].
In vivo	Brexpiprazole is able to ameliorate PCP- 191 induced cognitive deficits in mice, via 5-HT _{1A} receptors[2].
Cell Research	2.5 ng/ml of NGF(nerve growth factor) is used to study the potentiating effects of brexpiprazole on neurite outgrowth. Twenty-four hours after plating, the medium is replaced with DMEM medium containing 0.5% FBS and 1% penicillin-streptomycin with NGF (2.5 ng/ml), with or without brexpiprazole (0.001, 0.01, 0.1 or 1.0 μM). Four days after incubation with NGF (2.5 ng/ml) with or without drugs, morphometric analysis is performed on digitized images of live cells taken under phase-contrast illumination, with an inverted microscope linked to a camera. (Only for Reference)

Solubility Information

Solubility	H ₂ O: < 1 mg/mL (insoluble or slightly soluble), DMSO: 21 mg/mL (48.44 mM), Sonication is recommended. Ethanol: < 1 mg/mL (insoluble or slightly soluble), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3064 mL	11.5322 mL	23.0643 mL
5 mM	0.4613 mL	2.3064 mL	4.6129 mL
10 mM	0.2306 mL	1.1532 mL	2.3064 mL
50 mM	0.0461 mL	0.2306 mL	0.4613 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Ishima T, et al. Eur Neuropsychopharmacol. 2015, 25(4):505-511.

Yoshimi N, et al. Pharmacol Biochem Behav. 2014, 124:245-249.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481