


## Sendanolactone

## Chemical Properties

CAS No. : 64929-59-5  
 Formula: C<sub>30</sub>H<sub>42</sub>O<sub>4</sub>  
 Molecular Weight: 466.65  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sendanolactone possesses moderate cytotoxic activity against KB cell lines.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1429 mL	10.7147 mL	21.4293 mL
5 mM	0.4286 mL	2.1429 mL	4.2859 mL
10 mM	0.2143 mL	1.0715 mL	2.1429 mL
50 mM	0.0429 mL	0.2143 mL	0.4286 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481