


JTE 151A

Chemical Properties

CAS No. : 1608139-41-8
 Formula: C₂₆H₃₀ClN₃O₅
 Molecular Weight: 499.99
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	JTE 151A is a novel ROR γ inverse agonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.000 mL	10.0002 mL	20.0004 mL
5 mM	0.400 mL	2.000 mL	4.0001 mL
10 mM	0.200 mL	1.000 mL	2.000 mL
50 mM	0.040 mL	0.200 mL	0.400 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481