

## Methyl-PEG3-Ald

## Chemical Properties

CAS No. :	141981-62-6
Formula:	C <sub>9</sub> H <sub>18</sub> O <sub>5</sub>
Molecular Weight:	206.238
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Methyl-PEG3-Ald is a PEG-based linker that connects two essential ligands in PROTAC molecules, facilitating selective protein degradation via the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands connected by a linker: one ligand targets an E3 ubiquitin ligase, and the other targets the desired protein. They leverage the ubiquitin-proteasome system to achieve selective degradation of the target proteins[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8487 mL	24.2436 mL	48.4872 mL
5 mM	0.9697 mL	4.8487 mL	9.6974 mL
10 mM	0.4849 mL	2.4244 mL	4.8487 mL
50 mM	0.097 mL	0.4849 mL	0.9697 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481