


Tubulin inhibitor 23

Chemical Properties

CAS No. :	170488-57-0
Formula:	C ₂₆ H ₂₃ N ₀ O ₆ S
Molecular Weight:	477.53
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Tubulin inhibitor 23 is a potent tubulin inhibitor (IC ₅₀ : 4.8 μM) that exhibits antiangiogenic effects in a dose-dependent manner and induces apoptosis. This compound shows potential for leukemia research.
Targets(IC ₅₀)	Apoptosis,Others,Microtubule Associated

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0941 mL	10.4705 mL	20.9411 mL
5 mM	0.4188 mL	2.0941 mL	4.1882 mL
10 mM	0.2094 mL	1.0471 mL	2.0941 mL
50 mM	0.0419 mL	0.2094 mL	0.4188 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481