

## N1,N3-Bis(cyanomethyl)pseudouridine

## Chemical Properties

CAS No. :


Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	N1,N3-Bis(cyanomethyl)pseudouridine is a nucleoside derivative categorized as a C-nucleoside and an N-alkylated nucleoside.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481