

Cy3.5 diacid , SE

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3.5 diacid, SE is commonly used in biomolecular labeling, fluorescence imaging and other fluorescence biological analysis. Cy3.5 dyes and their derivatives are widely used to label polypeptides, proteins and nucleic acids.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481