


Rosmarinyl glucoside

Chemical Properties

CAS No. : 910028-78-3
 Formula: C₂₄H₂₆O₁₃
 Molecular Weight: 522.46
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description Rosmarinyl glucoside is a natural product for research related to life sciences. The catalog number is TN6523 and the CAS number is 910028-78-3.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.914 mL	9.5701 mL	19.1402 mL
5 mM	0.3828 mL	1.914 mL	3.828 mL
10 mM	0.1914 mL	0.957 mL	1.914 mL
50 mM	0.0383 mL	0.1914 mL	0.3828 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481