

5-endo-BCN-pentanoic acid

Chemical Properties

CAS No. : 2364591-80-8


Formula: C₁₆H₂₃NO₄

Molecular Weight: 293.363

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-endo-BCN-pentanoic acid, an alkyl chain-based PROTAC linker, facilitates the synthesis of PROTACs[1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs are molecules that consist of two distinct ligands connected by a linker: one ligand targets an E3 ubiquitin ligase, while the other binds to the target protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4088 mL	17.0439 mL	34.0878 mL
5 mM	0.6818 mL	3.4088 mL	6.8176 mL
10 mM	0.3409 mL	1.7044 mL	3.4088 mL
50 mM	0.0682 mL	0.3409 mL	0.6818 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481