


(Rac)-3'-Hydroxy simvastatin

Chemical Properties

CAS No. :	126313-98-2
Formula:	C ₂₅ H ₃₈ O ₆
Molecular Weight:	434.57
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(Rac)-3'-Hydroxy Simvastatin, a metabolite of Simvastatin, acts as a competitive inhibitor of HMG-CoA reductase, demonstrating a K _i value of 0.2 nM.
Targets(IC50)	Others, Drug Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3011 mL	11.5056 mL	23.0113 mL
5 mM	0.4602 mL	2.3011 mL	4.6023 mL
10 mM	0.2301 mL	1.1506 mL	2.3011 mL
50 mM	0.046 mL	0.2301 mL	0.4602 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481