


Eltrombopag methyl ester

Chemical Properties

CAS No. :	1246929-01-0
Formula:	C ₂₆ H ₂₄ N ₄ O ₄
Molecular Weight:	456.49
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Eltrombopag Methyl Ester, a derivative of Eltrombopag, functions as a thrombopoietin (Tpo) receptor agonist to enhance platelet production and is employed in the study of certain thrombocytopenia variants.
Targets(IC50)	Others,Thrombopoietin Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1906 mL	10.9531 mL	21.9063 mL
5 mM	0.4381 mL	2.1906 mL	4.3813 mL
10 mM	0.2191 mL	1.0953 mL	2.1906 mL
50 mM	0.0438 mL	0.2191 mL	0.4381 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481