


## Flurandrenolone Acetate

## Chemical Properties

CAS No. : 2802-11-1  
 Formula: C<sub>26</sub>H<sub>35</sub>F<sub>07</sub>  
 Molecular Weight: 478.55  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Flurandrenolone Acetate, a synthetic glucocorticoid steroid derivative of Flurandrenolide, is utilized in researching skin disorders including eczema and psoriasis [1].
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0896 mL	10.4482 mL	20.8965 mL
5 mM	0.4179 mL	2.0896 mL	4.1793 mL
10 mM	0.209 mL	1.0448 mL	2.0896 mL
50 mM	0.0418 mL	0.209 mL	0.4179 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481