


## Antiarrhythmic agent-1

## Chemical Properties

CAS No. : 136081-07-7  
 Formula: C<sub>25</sub>H<sub>27</sub>N<sub>3</sub>O<sub>4</sub>S  
 Molecular Weight: 465.56  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Antiarrhythmic agent-1 (example I), functioning as an antiarrhythmic agent and an IKr potassium channel blocker (IC <sub>50</sub> <1 μM), effectively inhibits cardiac arrhythmias.
Targets(IC <sub>50</sub> )	Others,Potassium Channel

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.148 mL	10.7398 mL	21.4795 mL
5 mM	0.4296 mL	2.148 mL	4.2959 mL
10 mM	0.2148 mL	1.074 mL	2.148 mL
50 mM	0.043 mL	0.2148 mL	0.4296 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481