

Lagociclovir valactate

Chemical Properties


CAS No. : 1001670-19-4

Formula: C₁₈H₂₅FN₆O₆

Molecular Weight: 440.43

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Lagociclovir valactate, a prodrug of Lagociclovir, is an orally active anti-HBV agent [1].
Targets(IC50)	Others,HBV

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2705 mL	11.3525 mL	22.7051 mL
5 mM	0.4541 mL	2.2705 mL	4.541 mL
10 mM	0.2271 mL	1.1353 mL	2.2705 mL
50 mM	0.0454 mL	0.2271 mL	0.4541 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481