

Siegesbeckialide I

Chemical Properties


CAS No. :

Formula: C₂₀H₂₈O₆

Molecular Weight: 364.43

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Siegesbeckialide I effectively suppresses LPS-induced NO production in RAW264.7 murine macrophages by directly interacting with IKK α / β .
Targets(IC50)	Others,I κ B/IKK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.744 mL	13.7201 mL	27.4401 mL
5 mM	0.5488 mL	2.744 mL	5.488 mL
10 mM	0.2744 mL	1.372 mL	2.744 mL
50 mM	0.0549 mL	0.2744 mL	0.5488 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481