


## Dihydroaltenuene B

## Chemical Properties

CAS No. : 887751-89-5  
 Formula: C<sub>15</sub>H<sub>18</sub>O<sub>6</sub>  
 Molecular Weight: 294.3  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Dihydroaltenuene B acts as an effective mushroom tyrosinase inhibitor, exhibiting an inhibition concentration (IC <sub>50</sub> ) of 38.33 μM. It demonstrates hydrogen bonding interactions with the His244, Met280, and Gly281 residues of tyrosinase through its 3-OH and 4'-OH groups.
Targets(IC <sub>50</sub> )	Others,Antifungal,Tyrosinase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3979 mL	16.9895 mL	33.9789 mL
5 mM	0.6796 mL	3.3979 mL	6.7958 mL
10 mM	0.3398 mL	1.6989 mL	3.3979 mL
50 mM	0.068 mL	0.3398 mL	0.6796 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481