


α,β -Methylene-ATP

Chemical Properties

CAS No. :	7292-42-4
Formula:	C ₁₁ H ₁₈ N ₅ O ₁₂ P ₃
Molecular Weight:	505.21
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	α,β -Methylene ATP, a phosphonic analog of ATP, functions as a ligand for P2X ₃ and P2X ₇ receptors, and acts primarily as a selective agonist for P2X ₁ and P2X ₃ receptors, exhibiting minimal to no activity on P2X ₂ , P2X ₄ through P2X ₇ receptors.
Targets(IC50)	Others,P2X Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9794 mL	9.8969 mL	19.7937 mL
5 mM	0.3959 mL	1.9794 mL	3.9587 mL
10 mM	0.1979 mL	0.9897 mL	1.9794 mL
50 mM	0.0396 mL	0.1979 mL	0.3959 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481