


## IL-17 modulator 5

## Chemical Properties

CAS No. : 2724206-27-1  
 Formula: C<sub>28</sub>H<sub>23</sub>F<sub>6</sub>N<sub>9</sub>O<sub>2</sub>  
 Molecular Weight: 631.53  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	IL-17 modulator 5 is a IL-17 inhibitor, with an IC <sub>50</sub> of 1 nM .
Targets(IC <sub>50</sub> )	Others,Interleukin

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5835 mL	7.9173 mL	15.8346 mL
5 mM	0.3167 mL	1.5835 mL	3.1669 mL
10 mM	0.1583 mL	0.7917 mL	1.5835 mL
50 mM	0.0317 mL	0.1583 mL	0.3167 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481