


AT1R antagonist 1

Chemical Properties

CAS No. :	2848719-80-0
Formula:	C ₂₈ H ₃₂ F ₃ N ₅ O ₄ S ₂
Molecular Weight:	623.71
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	AT1R Antagonist 1 (compound 10) serves as a potent, selective ligand for AT1R with considerable affinity, demonstrated by a K _i value of 8.5 nM [1].
Targets(IC50)	RAAS,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6033 mL	8.0165 mL	16.0331 mL
5 mM	0.3207 mL	1.6033 mL	3.2066 mL
10 mM	0.1603 mL	0.8017 mL	1.6033 mL
50 mM	0.0321 mL	0.1603 mL	0.3207 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481