


## GLP-1R agonist 17

## Chemical Properties

CAS No. :	2749609-28-5
Formula:	C <sub>28</sub> H <sub>26</sub> ClFN <sub>4</sub> O <sub>4</sub> S
Molecular Weight:	569.05
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	GLP-1R agonist 17 is a potent GLP-1 receptor agonist with excellent efficacy in activating the GLP-1 receptor, primarily used in cardiovascular metabolic disease research.
Targets(IC50)	Others, Glucagon Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7573 mL	8.7866 mL	17.5731 mL
5 mM	0.3515 mL	1.7573 mL	3.5146 mL
10 mM	0.1757 mL	0.8787 mL	1.7573 mL
50 mM	0.0351 mL	0.1757 mL	0.3515 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481