


## Loxoprofenol-SRS

## Chemical Properties

CAS No. :	83648-76-4
Formula:	C <sub>15</sub> H <sub>20</sub> O <sub>3</sub>
Molecular Weight:	248.32
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Loxoprofenol-SRS, a new intravenous non-steroidal anti-inflammatory drug (NSAID) and the active metabolite of Loxoprofen, demonstrates significantly enhanced anti-inflammatory and analgesic (pain-relieving) effects.
Targets(IC50)	Others, Drug Metabolite

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0271 mL	20.1353 mL	40.2706 mL
5 mM	0.8054 mL	4.0271 mL	8.0541 mL
10 mM	0.4027 mL	2.0135 mL	4.0271 mL
50 mM	0.0805 mL	0.4027 mL	0.8054 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481