

## Maydispenoid B

## Chemical Properties

CAS No. :	3024647-69-3
Formula:	C <sub>25</sub> H <sub>34</sub> O <sub>4</sub>
Molecular Weight:	398.54
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Maydispenoid B, a potent immunosuppressor, effectively inhibits the proliferation of murine splenocytes activated by anti-CD3/anti-CD28 monoclonal antibodies (mAbs) and lipopolysaccharide.
Targets(IC50)	Others, Immunology/Inflammation related, Integrin

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5092 mL	12.5458 mL	25.0916 mL
5 mM	0.5018 mL	2.5092 mL	5.0183 mL
10 mM	0.2509 mL	1.2546 mL	2.5092 mL
50 mM	0.0502 mL	0.2509 mL	0.5018 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481