

Dynamamin IN-1

Chemical Properties


CAS No. : 1345853-50-0

Formula: C₂₃H₂₄N₂O

Molecular Weight: 344.45

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dynamamin IN-1, a potent dynamamin inhibitor, has an IC ₅₀ value of 1.0 μM.
Targets(IC50)	Others,Dynamamin

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9032 mL	14.5159 mL	29.0318 mL
5 mM	0.5806 mL	2.9032 mL	5.8064 mL
10 mM	0.2903 mL	1.4516 mL	2.9032 mL
50 mM	0.0581 mL	0.2903 mL	0.5806 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481