

## Sycosterol A

## Chemical Properties


CAS No. :

Formula: C27H48O9S

Molecular Weight: 548.73

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Sycosterol A is an $\alpha$ -synuclein ( $\alpha$ -syn) inhibitory sterol.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8224 mL	9.1119 mL	18.2239 mL
5 mM	0.3645 mL	1.8224 mL	3.6448 mL
10 mM	0.1822 mL	0.9112 mL	1.8224 mL
50 mM	0.0364 mL	0.1822 mL	0.3645 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481