


## Dictyostatin

## Chemical Properties

CAS No. :	156312-07-1
Formula:	C <sub>32</sub> H <sub>52</sub> O <sub>6</sub>
Molecular Weight:	532.75
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Dictyostatin ((-)-Dictyostatin; Dictyostatin 1) is a potent microtubule-stabilizing and anti-cancer agent exhibiting antiproliferative activity, with potential applications in tauopathy research.
Targets(IC50)	Others, Microtubule Associated

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8771 mL	9.3853 mL	18.7705 mL
5 mM	0.3754 mL	1.8771 mL	3.7541 mL
10 mM	0.1877 mL	0.9385 mL	1.8771 mL
50 mM	0.0375 mL	0.1877 mL	0.3754 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481