


(S)-JDQ-443

## Chemical Properties

CAS No. : 2653994-10-4  
 Formula: C<sub>29</sub>H<sub>28</sub>ClN<sub>7</sub>O  
 Molecular Weight: 526.03  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	(S)-JDQ-443 is an isomer of JDQ-443, which is an orally active, potent, selective, and covalent KRAS G12C inhibitor, showing antitumor activity [1] [2].
Targets(IC50)	Ras,Kras

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.901 mL	9.5052 mL	19.0103 mL
5 mM	0.3802 mL	1.901 mL	3.8021 mL
10 mM	0.1901 mL	0.9505 mL	1.901 mL
50 mM	0.038 mL	0.1901 mL	0.3802 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481