

Aflibercept

Chemical Properties

CAS No. : 862111-32-8

Formula:

Molecular Weight:

Store at low temperature

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Aflibercept is a recombinant fusion protein antibody composed of the Ig domains of VEGFR1 and VEGFR2 and the Fc domain of human IgG1. Aflibercept inhibits VEGF signaling and is a commonly used tool in biological research for suppressing pathological angiogenesis.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481