


STAT3-IN-14

Chemical Properties

CAS No. :	123297-90-5
Formula:	C ₁₄ H ₁₀ O ₅
Molecular Weight:	258.23
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	STAT3-IN-14 is a STAT3 inhibitor and has STAT3 phosphorylation inhibitory activity. STAT3-IN-14 can directly bind to the hinge region of STAT3 .
Targets(IC50)	Others,STAT

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8725 mL	19.3626 mL	38.7252 mL
5 mM	0.7745 mL	3.8725 mL	7.745 mL
10 mM	0.3873 mL	1.9363 mL	3.8725 mL
50 mM	0.0775 mL	0.3873 mL	0.7745 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481