


Laburnetin

Chemical Properties

CAS No. :	166375-17-3
Formula:	C ₂₀ H ₁₈ O ₆
Molecular Weight:	354.35
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Laburnetin, an isoflavone antibacterial agent, exhibits antibacterial activity against fungi and <i>S. vesicarium</i> , and enhances the susceptibility of Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA) strains to Methicillin. It is also effective in controlling pests of cultivated species [1] [2] [3].
Targets(IC50)	Others,Antibacterial

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8221 mL	14.1103 mL	28.2207 mL
5 mM	0.5644 mL	2.8221 mL	5.6441 mL
10 mM	0.2822 mL	1.411 mL	2.8221 mL
50 mM	0.0564 mL	0.2822 mL	0.5644 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481