


Gepefrine

Chemical Properties

CAS No. :	18840-47-6
Formula:	C ₉ H ₁₃ NO
Molecular Weight:	151.21
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Gepefrine, an orally active pressor and sympathomimetic agent, enhances the initial orthostatic regulation of arterial pressure.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.6133 mL	33.0666 mL	66.1332 mL
5 mM	1.3227 mL	6.6133 mL	13.2266 mL
10 mM	0.6613 mL	3.3067 mL	6.6133 mL
50 mM	0.1323 mL	0.6613 mL	1.3227 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481