


BTK-IN-22

Chemical Properties

CAS No. : 2573048-10-7
 Formula: C₂₆H₂₆N₆O₂
 Molecular Weight: 454.52
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	BTK-IN-22, a BTK inhibitor (IC 50: 0.9 nM), also exhibits inhibition against BLX and BMX with IC 50 values of 1.4 and 1.2 nM, respectively. Demonstrating enhanced kinase selectivity, BTK-IN-22 outperforms Ibrutinib in this regard [1].
Targets(IC50)	Others,BTK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2001 mL	11.0006 mL	22.0012 mL
5 mM	0.440 mL	2.2001 mL	4.4002 mL
10 mM	0.220 mL	1.1001 mL	2.2001 mL
50 mM	0.044 mL	0.220 mL	0.440 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481