


Vonifimod

Chemical Properties

CAS No. :	2118903-18-5
Formula:	C17H34N4O2
Molecular Weight:	326.48
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Vonifimod is an immunomodulator .
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.063 mL	15.3149 mL	30.6297 mL
5 mM	0.6126 mL	3.063 mL	6.1259 mL
10 mM	0.3063 mL	1.5315 mL	3.063 mL
50 mM	0.0613 mL	0.3063 mL	0.6126 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481