


AP 24149

Chemical Properties

CAS No. : 926922-29-4
 Formula: C₂₃H₂₄N₅O_P
 Molecular Weight: 417.44
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	AP 24149, a potent dual inhibitor targeting Src and Abl, exhibits IC ₅₀ values of 9.1 nM for Src and 3.6 nM for Abl, respectively.
Targets(IC ₅₀)	Others,Bcr-Abl,Src

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3956 mL	11.9778 mL	23.9555 mL
5 mM	0.4791 mL	2.3956 mL	4.7911 mL
10 mM	0.2396 mL	1.1978 mL	2.3956 mL
50 mM	0.0479 mL	0.2396 mL	0.4791 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481