

BI-2540

## Chemical Properties


CAS No. : 875145-22-5

Formula: C<sub>24</sub>H<sub>15</sub>ClF<sub>5</sub>NO<sub>5</sub>

Molecular Weight: 527.83

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	BI-2540 is a non-nucleoside reverse transcriptase (NNRT) inhibitor for HIV.
Targets(IC50)	Others,HIV Protease

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8945 mL	9.4727 mL	18.9455 mL
5 mM	0.3789 mL	1.8945 mL	3.7891 mL
10 mM	0.1895 mL	0.9473 mL	1.8945 mL
50 mM	0.0379 mL	0.1895 mL	0.3789 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481