

PG 34

Chemical Properties


CAS No. : 1681056-62-1

Formula: C₂₈H₂₅N₃O₅S₂

Molecular Weight: 547.65

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	PG 34 is a structural analogue of OSMI-1 and a weak inhibitor of O-GlcNAc translocase (OGT) with an IC ₅₀ value of 68 μM.
Targets(IC ₅₀)	Others,Transferase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.826 mL	9.1299 mL	18.2598 mL
5 mM	0.3652 mL	1.826 mL	3.652 mL
10 mM	0.1826 mL	0.913 mL	1.826 mL
50 mM	0.0365 mL	0.1826 mL	0.3652 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481