


## Compound EN300-219925

## Chemical Properties

CAS No. : 1803595-30-3  
 Formula: C<sub>14</sub>H<sub>12</sub>FNO<sub>2</sub>S  
 Molecular Weight: 277.31  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound EN300-219925 is a high-purity biochemical reagent with CAS number 1803595-30-3, suitable for biochemical experiments and drug synthesis research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6061 mL	18.0304 mL	36.0607 mL
5 mM	0.7212 mL	3.6061 mL	7.2121 mL
10 mM	0.3606 mL	1.803 mL	3.6061 mL
50 mM	0.0721 mL	0.3606 mL	0.7212 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481