


3-Isomangostin hydrate

Chemical Properties

CAS No. :	26063-96-7
Formula:	C ₂₄ H ₂₈ O ₇
Molecular Weight:	428.47
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-Isomangostin hydrate is a natural product for research related to life sciences. The catalog number is TN2954 and the CAS number is 26063-96-7.
Targets(IC50)	Others,Antibacterial

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3339 mL	11.6694 mL	23.3389 mL
5 mM	0.4668 mL	2.3339 mL	4.6678 mL
10 mM	0.2334 mL	1.1669 mL	2.3339 mL
50 mM	0.0467 mL	0.2334 mL	0.4668 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481