


## Agastenol

## Chemical Properties

CAS No. :	413615-34-6
Formula:	C <sub>27</sub> H <sub>26</sub> O <sub>8</sub>
Molecular Weight:	478.497
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Agastenol is a natural product that can be used as a reference standard. The CAS number of Agastenol is 413615-34-6.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0899 mL	10.4493 mL	20.8986 mL
5 mM	0.418 mL	2.0899 mL	4.1797 mL
10 mM	0.209 mL	1.0449 mL	2.0899 mL
50 mM	0.0418 mL	0.209 mL	0.418 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481