

3-Epiglochidiol diacetate

Chemical Properties


CAS No. : 6587-37-7

Formula: C₃₄H₅₄O₄

Molecular Weight: 526.79

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-Epiglochidiol diacetate is a natural product for research related to life sciences. The catalog number is TN2930 and the CAS number is 6587-37-7.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8983 mL	9.4914 mL	18.9829 mL
5 mM	0.3797 mL	1.8983 mL	3.7966 mL
10 mM	0.1898 mL	0.9491 mL	1.8983 mL
50 mM	0.038 mL	0.1898 mL	0.3797 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481