


DK-139

Chemical Properties

CAS No. : 1426059-23-5
 Formula: C₁₈H₁₈O₅
 Molecular Weight: 314.33
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	DK-139 is an inhibitor of NF-kappaB-mediated GROalpha expression, suppressing the TNFalpha-induced invasive capability of MDA-MB-231 human breast cancer cells.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1814 mL	15.9068 mL	31.8137 mL
5 mM	0.6363 mL	3.1814 mL	6.3627 mL
10 mM	0.3181 mL	1.5907 mL	3.1814 mL
50 mM	0.0636 mL	0.3181 mL	0.6363 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481