

Valtropine

Chemical Properties

CAS No. :	495-82-9
Formula:	C13H23NO2
Molecular Weight:	225.33
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Valtropine is a useful organic compound for research related to life sciences and the catalog number is T124775.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4379 mL	22.1897 mL	44.3794 mL
5 mM	0.8876 mL	4.4379 mL	8.8759 mL
10 mM	0.4438 mL	2.219 mL	4.4379 mL
50 mM	0.0888 mL	0.4438 mL	0.8876 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481