


Asp-Asp-Asp-Asp-Asp

Chemical Properties

CAS No. :	124219-00-7
Formula:	C ₂₀ H ₂₇ N ₅ O ₁₆
Molecular Weight:	593.45
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Asp-Asp-Asp-Asp-Asp is a peptide consisting of five Asp (aspartic acid) residues.
-------------	---

Solubility Information

Solubility	H ₂ O: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6851 mL	8.4253 mL	16.8506 mL
5 mM	0.337 mL	1.6851 mL	3.3701 mL
10 mM	0.1685 mL	0.8425 mL	1.6851 mL
50 mM	0.0337 mL	0.1685 mL	0.337 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481