


## (2R)-1-Triphenylmethyl-2-pyrrolidinecarboxylic acid methyl ester

### Chemical Properties

CAS No. :	1033699-77-2
Formula:	C <sub>25</sub> H <sub>25</sub> NO <sub>2</sub>
Molecular Weight:	371.47
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



### Biological Description

Description	(2R)-1-Triphenylmethyl-2-pyrrolidinecarboxylic acid methyl ester belongs to Intermediates and Building Blocks - Amino acid.
Targets(IC50)	Others

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.692 mL	13.460 mL	26.9201 mL
5 mM	0.5384 mL	2.692 mL	5.384 mL
10 mM	0.2692 mL	1.346 mL	2.692 mL
50 mM	0.0538 mL	0.2692 mL	0.5384 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481