


L-Azidohomoalanine hydrochloride

Chemical Properties

CAS No. :	942518-29-8
Formula:	C ₄ H ₉ ClN ₄ O ₂
Molecular Weight:	180.59
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	L-Azidohomoalanine hydrochloride (H-L-Aha-OH hydrochloride) is an alkyl chain-based PROTAC linker.
Targets(IC50)	PROTAC Linker

Solubility Information

Solubility	H ₂ O: 250 mg/mL (1384.35 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.5374 mL	27.687 mL	55.3741 mL
5 mM	1.1075 mL	5.5374 mL	11.0748 mL
10 mM	0.5537 mL	2.7687 mL	5.5374 mL
50 mM	0.1107 mL	0.5537 mL	1.1075 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481