

PDL-1 cpd 10

Chemical Properties


CAS No. : 2767424-13-3

Formula: C₂₁H₂₃N₅O₂

Molecular Weight: 377.44

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	PDL-1 cpd 10 is a novel small molecule PD-L1 inhibitor.
Targets(IC50)	Apoptosis,PD-1/PD-L1

Solubility Information

Solubility	DMSO: 55 mg/mL (145.72 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (5.3 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6494 mL	13.2471 mL	26.4943 mL
5 mM	0.5299 mL	2.6494 mL	5.2989 mL
10 mM	0.2649 mL	1.3247 mL	2.6494 mL
50 mM	0.053 mL	0.2649 mL	0.5299 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Russomanno P, et al. Interfering with the Tumor-Immune Interface: Making Way for Triazine-Based Small Molecules as Novel PD-L1 Inhibitors. J Med Chem. 2021 Nov 11;64(21):16020-16045.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481