

Biotin-PEG-Biotin (MW 1000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Biotin-PEG-Biotin (MW 1000) is a PEG-derived PROTAC linker utilized in PROTAC synthesis [1]. |
| Targets(IC50) | Others,PROTAC Linker |
| In vitro | PROTACs consist of two distinct ligands linked together: one binds to an E3 ubiquitin ligase, and the other targets a specific protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins [1]. |

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481