


## Shizukanolide H

## Chemical Properties

CAS No. : 1136932-34-7  
 Formula: C<sub>17</sub>H<sub>20</sub>O<sub>5</sub>  
 Molecular Weight: 304.34  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Shizukanolide H is a natural product of Chloranthus, Chloranthaceae. The catalog number is TN5009 and the CAS number is 1136932-34-7. Shizukanolide H can be used as a reference standard.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2858 mL	16.429 mL	32.858 mL
5 mM	0.6572 mL	3.2858 mL	6.5716 mL
10 mM	0.3286 mL	1.6429 mL	3.2858 mL
50 mM	0.0657 mL	0.3286 mL	0.6572 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481