


PCSK9 modulator-4

Chemical Properties

CAS No. :	2476490-20-5
Formula:	C ₁₇ H ₁₁ F ₂ N ₃ O
Molecular Weight:	311.29
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	PCSK9 modulator-4 (Compound 21) is a potent PCSK9 modulator (EC ₅₀ = 0.15 nM) utilized in hyperlipidemia research, targeting the recently validated low-density lipoprotein cholesterol (LDL-C) lowering mechanism [1].
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2124 mL	16.0622 mL	32.1244 mL
5 mM	0.6425 mL	3.2124 mL	6.4249 mL
10 mM	0.3212 mL	1.6062 mL	3.2124 mL
50 mM	0.0642 mL	0.3212 mL	0.6425 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481